


Best™ Cyclotron Systems


TEAMBEST GLOBAL®



www.teambest.com




Best BG-95 Sub-Compact Self-Shielded Cyclotron with Optional Second Chemistry Module




Best Particle Therapy 400 MeV ion Rapid Cycling Medical Synchrotron (iRCMS)


NEW! Best Model 200p Cyclotron—Variable Energy Proton Beam for Radiation Therapy*



*Patent Pending




Best Model B35ADP Alpha/Deuteron/Proton Cyclotron




Best Model 70p Cyclotron for Radioisotopes Production and Research

Best Model 6-15 MeV Compact High Current/Variable Energy Proton Cyclotron



X-Beam™ Robotic Radiosurgery System



X-Beam™ Image-Guided Multi-Energy Linac System

NEW! B-11	11 MeV	• For radioisotopes production and research
B100 Cyclotron	7.5 MeV	• Capable of producing: ^{18}F and Na^{18}F • Single or batch dose production • Integrated self-shielded cyclotron, chemistry module and FDG QC module • Complete production lab in a 5 x 5 meter area
BG-95 Cyclotron	1-9.5 MeV	• Low energy, self-shielded compact system capable of producing: ^{18}F , Na^{18}F , ^{18}F -MISO, ^{18}F -FLT, ^{18}F -Choline, ^{18}F -DOPA, ^{18}F -PSMA, ^{13}N and ^{68}Ga
Best Cyclotrons	1-3 MeV	• Deuterons for materials analysis*
	70-200 MeV	• For Proton Therapy*
	3-90 MeV	• High current proton beams for neutron production and delivery*
B6-15 Cyclotron	1-15 MeV	• Proton only, capable of high current up to 1000 Micro Amps, for medical radioisotopes
B25 Cyclotron	20, 15-25 MeV	• Proton only, capable of high current up to 1000 Micro Amps, for medical radioisotopes
B25u-35adp Cyclotron	25-35 MeV	• Proton or alpha/deuteron/proton, capable of high current up to 1000 Micro Amps, for medical radioisotopes
B35 Cyclotron	35 MeV	• Proton only system for medical radioisotopes production
B70/70adp Cyclotron	35-70 MeV	• Proton only or alpha/deuteron/proton systems, capable of high current up to 1000 Micro Amps, for medical radioisotopes

*Some products are under development and not available for sale currently.

FLASH THERAPY





E-Beam™ Robotic IORT Linac System

TeamBest Global Companies ©2024

Best Medical International, 7643 Fullerton Road, Springfield, VA 22153
tel: 703 451 2378 800 336 4970 www.besttotalsolutions.com www.teambest.com

Best Theratronics Ltd., 413 March Road, Ottawa, Ontario K2K 0E4 Canada
tel: 613 591 2100 866 792 8598 www.theratronics.com



Best GammaBeamTM 100/300 EquinoxTM Teletherapy System with AvanzaTM 6D Patient Positioning Table

With **NEW** Multi-Leaf Collimator for 80 and 100 cm SAD units—
IMRT, IGRT, SRS, SBRT and Tomotherapy capable with ActiveRx

UPGRADE KIT for all old Theratron units, 80 or 100 cm including IMRT capabilities w/built in or external MLC



UPGRADE includes:

- Removing all old controls, electronics and installing a new control system and covers
- Replacing the old collimator system with the new Equinox collimator
- Replacing the old treatment table with the new AvanzaTM Table
- Retaining the head rotation capability is optional

UPGRADE features:

- Calculated Arc Speed
- Graphical Control System
- Asymmetric Jaws (optional)
- R&V System Ready (optional)
- Service Log Files
- On-Board Verification
- Motorized Wedge (optional)
- Collision Detection (optional)



BestTM Gammacell® 1000/3000
Blood & Research Irradiator



BestTM Raycell
Mk2 Blood
& Research
Irradiator



BestTM Raycell Mk1 Blood
& Research Irradiator



Best[™] medical international

A TEAMBEST GLOBAL COMPANY


Best[™] Integrated Brachytherapy Solutions

You don't go to a shoe store and cut your feet to fit the shoe they have in stock ... then why compromise when it comes to improving the clinical outcome for your patients?



Best[™] Universal Chair/Table

Best[™] HDR Remote Afterloader



BestMick Applicator

BestMick 10/15/20/25 Cartridges

Best[™] Kobold Tandem & Ring Applicator


Best[™] Kobold Fletcher Tandem & Ovoid Applicator

Best[™] Kobold Henschke Tandem & Ovoid Applicator

Best[™] Compact Sonalis[™] Ultrasound Imaging System



Best[™] X-ray Marker Cables




mobileMOSFET Wireless Dosimetry System



Portable Dosimetry System

Best[™] Templates



Best[™] Stepper/Stabilizer



Best[™] Treatment Planning System



Best[™] Esophageal Brachytherapy Applicator




Best[™] Breast Double-Balloon Brachytherapy Applicator



Best[™] Localization Needles with I-125 Seeds




Best[™] Fiducial Markers



Best[™] Iodine-125 Seed

Best[™] Palladium-103 Seed

Best[™] Flexi Needles



15G Flexi Needle

16G Flexi Needle

17G Flexi Needle

15G Flexi Needle with SS Sleeve (Click Fit)

Flexi Needle with Female Hub

Huestis CVO-2000 Warming Oven for Thermoplastics



We ship your order within 24 hours, 7 days a week, anywhere in the world!
Standard or customized marker configuration made to your exact specifications, sterile or non-sterile!

Best Medical International, 7643 Fullerton Road, Springfield, VA 22153

tel: 703 451 2378 800 336 4970 www.besttotalsolutions.com www.teambest.com

*Certain products shown are not available for sale currently. TeamBest Global Companies ©2024

Best[™] nomos[®] Compact Sonalis[®] Ultrasound System



The Best[™] Sonalis[®] Ultrasound Imaging System provides superior visualization of HDR, LDR, RF or Cryosurgical procedures. Our patented SimulView[™] Technology provides simultaneous “live” views of the prostate in both planes, thereby increasing treatment accuracy and precision.

Additional Features

- Capable of supporting over 20 different probe configurations, this system can be used in almost every discipline
- 21.5" touch screen technology
- Superior HD image resolution for improved procedure accuracy, speed and physician confidence level
- Patented probe design
- Sagittal array provides for 140 mm length of view encompassing the bladder, prostate and perineal tissue
- Advanced modular software design provides for future upgrade path
- Convertible to desktop or portable without losing any functionality
- Dual battery backup built-in to protect from power loss



*Certain products shown are not available for sale currently.

Sonalis[®] Transducers

We offer probes for the most common types of procedures and imaging, all with the enhanced image quality and high resolution that you would expect from the **BEST** name in ultrasounds!


- 8L2A Linear Array
- 12L5A Linear Array
- 14L3 Linear Array
- 15LW4 Linear Array
- 15LA Linear Array
- 15L4A Linear Array
- 16L5 Linear Array
- 8V3 Phased Array
- 4V2A Phased Array
- 5C2A Curved Array
- 9MC3 Curved Array
- 8EC4A Endocavity
- XY-BI-Plane Phased Array
- 10EC4A Endocavity
- 10BP4 Bi-Plane
- 8BP4 Bi-Plane
- 8TE3 Trans-esophageal
- Pedoff
- 16HL7 High Frequency Linear Array

TeamBest Global Companies ©2024

Best Nomos • 1 Best Drive • Pittsburgh, PA 15202 USA • tel 412 312 6700 or 800 70 NOMOS • fax 412 312 6701 • www.nomos.com

AccuBoost[®] DIGITAL PARTIAL BREAST BRACHYTHERAPY SOLUTION

AccuBoost[®] Digital—A digital mammography-guided Non-Invasive Breast Brachytherapy (NIBB) equipment solution. This complete, turn-key solution provides all the components needed for this advanced treatment, bundled together in a streamlined offering.



HUESTIS MEDICAL
A TEAMBEST GLOBAL COMPANY

One Industrial Way | Tyngsboro, MA 01879
1-800-972-9222

www.accuboot.com
www.huestis.com
www.teambest.com



Welcome to Our World

TeamBest Global

www.teambest.com

*Ready to serve all your
needs globally*

NON-PROFIT

- **Best Cure Foundation**
www.bestcure.md
- **Brachytherapy Research & Educational Foundation**
www.brachytherapy.org
- **Global Best Cure Foundation**
www.globalbestcure.org

FOR-PROFIT

- **Arplay Medical**
www.arplay.com
- **Best ABT**
www.bestabt.com
- **Best Automation & Robotics**
www.teambestautomation.com

- **Best Cyclotron Systems, Inc.**
www.bestcyclotron.com
- **Best Dosimetry Services**
www.bestdosimetry.com
- **Best Entertainment**
www.bestentertainment.ca
- **Best Medical Canada**
www.bestmedical.ca or
www.mosfet.ca
- **Best Medical Capital, Inc.**
www.teambest.com
- **Best Medical International**
www.bestmedical.com
- **Best NOMOS**
www.nomos.com

- **Best Particle Therapy**
www.bestproton.com
- **Best Theratronics**
www.theratronics.com
- **Best Vascular (Novoste)**
www.bestvascular.com or
www.novoste.com
- **CNMC Company**
www.cnmcco.com
- **Huestis Medical**
www.huestis.com
- **Kitsault Energy**
www.kitsaultenergy.com

TEAMBEST GLOBAL®



www.teambest.com

Best Medical International
7643 Fullerton Road, Springfield, VA 22153
tel: 703 451 2378 800 336 4970

TeamBest Global Companies ©2024

TEAMBEST GLOBAL®



www.teambest.com