


Best™ Cyclotron Systems


TEAMBEST GLOBAL®



www.teambest.com





Best BG-95 Sub-Compact Self-Shielded Cyclotron with Optional Second Chemistry Module




Best Model B35ADP Alpha/Deuteron/Proton Cyclotron

Best Model 6-15 MeV Compact High Current/Variable Energy Proton Cyclotron




Best Particle Therapy 400 MeV ion Rapid Cycling Medical Synchrotron (iRCMS)



Best Model 70p Cyclotron for Radioisotopes Production and Research




X-Beam™ Robotic Radiosurgery System



X-Beam™ Image-Guided Multi-Energy Linac System

NEW! Best Model 200p Cyclotron—Variable Energy Proton Beam for Radiation Therapy*



*Patent Pending

NEW! B-11	11 MeV	• For radioisotopes production and research
B100 Cyclotron	7.5 MeV	• Capable of producing: ¹⁸ FDG and Na ¹⁸ F • Single or batch dose production • Integrated self-shielded cyclotron, chemistry module and FDG QC module • Complete production lab in a 5 x 5 meter area
BG-95 Cyclotron	1-9.5 MeV	• Low energy, self-shielded compact system capable of producing: ¹⁸ FDG, Na ¹⁸ F, ¹⁸ F-MISO, ¹⁸ FLT, ¹⁸ F-Choline, ¹⁸ F-DOPA, ¹⁸ F-PSMA, ¹³ N and ⁶⁸ Ga
Best Cyclotrons	1-3 MeV	• Deuterons for materials analysis*
	70-200 MeV	• For Proton Therapy*
	3-90 MeV	• High current proton beams for neutron production and delivery*
B6-15 Cyclotron	1-15 MeV	• Proton only, capable of high current up to 1000 Micro Amps, for medical radioisotopes
B25 Cyclotron	20, 15-25 MeV	• Proton only, capable of high current up to 1000 Micro Amps, for medical radioisotopes
B25u-35adp Cyclotron	25-35 MeV	• Proton or alpha/deuteron/proton, capable of high current up to 1000 Micro Amps, for medical radioisotopes
B35 Cyclotron	35 MeV	• Proton only system for medical radioisotopes production
B70/70adp Cyclotron	35-70 MeV	• Proton only or alpha/deuteron/proton systems, capable of high current up to 1000 Micro Amps, for medical radioisotopes

*Some products are under development and not available for sale currently.



E-Beam™ Robotic IORT Linac System

TeamBest Global Companies ©2024

Best Medical International, 7643 Fullerton Road, Springfield, VA 22153
tel: 703 451 2378 800 336 4970 www.besttotalsolutions.com www.teambest.com

Best Theratronics Ltd., 413 March Road, Ottawa, Ontario K2K 0E4 Canada
tel: 613 591 2100 866 792 8598 www.theratronics.com



Best™ *medical international*

A TEAMBEST GLOBAL COMPANY

TEAMBEST GLOBAL®



www.teambest.com



Best GammaBeam™ 100/300 Equinox™ Teletherapy System with Avanza™ 6D Patient Positioning Table

With **NEW** Multi-Leaf Collimator for 80 and 100 cm SAD units—IMRT, IGRT, SRS, SBRT and Tomotherapy capable with ActiveRx



Best™ Gammacell® 1000/3000 Blood & Research Irradiator

UPGRADE KIT for all old Theratron units, 80 or 100 cm including IMRT capabilities w/built in or external MLC



UPGRADE includes:

- Removing all old controls, electronics and installing a new control system and covers
- Replacing the old collimator system with the new Equinox collimator
- Replacing the old treatment table with the new Avanza™ Table
- Retaining the head rotation capability is optional

UPGRADE features:

- Calculated Arc Speed
- Graphical Control System
- Asymmetric Jaws (optional)
- R&V System Ready (optional)
- Service Log Files
- On-Board Verification
- Motorized Wedge (optional)
- Collision Detection (optional)



Best™ Raycell Mk2 Blood & Research Irradiator



Best™ Raycell Mk1 Blood & Research Irradiator

Best Theratronics Ltd., 413 March Road, Ottawa, Ontario K2K 0E4 Canada
tel: 613 591 2100 866 792 8598 www.theratronics.com www.teambest.com

*Certain products shown are not available for sale currently. TeamBest Global Companies ©2024



Best™ medical international

A TEAMBEST GLOBAL COMPANY

Best™ Integrated Brachytherapy Solutions

You don't go to a shoe store and cut your feet to fit the shoe they have in stock ... then why compromise when it comes to improving the clinical outcome for your patients?



Best™ Universal Chair/Table


Best™ HDR Remote Afterloader



BestMick Applicator



BestMick 10/15/20/25 Cartridges



Best™ Kobold Tandem & Ring Applicator


Best™ Kobold Fletcher Tandem & Ovoid Applicator

Best™ Kobold Henschke Tandem & Ovoid Applicator

Best™ Compact Sonalis™ Ultrasound Imaging System



Best™ X-ray Marker Cables




mobileMOSFET Wireless Dosimetry System




Portable Dosimetry System

Best™ Templates




Best™ Stepper/Stabilizer




Best™ Treatment Planning System



Best™ Esophageal Brachytherapy Applicator




Best™ Breast Double-Balloon Brachytherapy Applicator



Best™ Localization Needles with I-125 Seeds

Best™ Fiducial Markers




Best™ Iodine-125 Seed



Best™ Palladium-103 Seed



Best™ Flexi Needles



15G Flexi Needle

16G Flexi Needle

17G Flexi Needle

15G Flexi Needle with SS Sleeve (Click Fit)

Flexi Needle with Female Hub



Huestis CVO-2000 Warming Oven for Thermoplastics



We ship your order within 24 hours, 7 days a week, anywhere in the world!

Standard or customized marker configuration made to your exact specifications, sterile or non-sterile!

Best Medical International, 7643 Fullerton Road, Springfield, VA 22153

tel: 703 451 2378 800 336 4970 www.besttotalsolutions.com www.teambest.com

*Certain products shown are not available for sale currently. TeamBest Global Companies ©2024

Best[™] nomos[®] Compact Sonalis[®] Ultrasound System



The Best[™] Sonalis[®] Ultrasound Imaging System provides superior visualization of HDR, LDR, RF or Cryosurgical procedures. Our patented SimuView[™] Technology provides simultaneous “live” views of the prostate in both planes, thereby increasing treatment accuracy and precision.



Additional Features

- Capable of supporting over 20 different probe configurations, this system can be used in almost every discipline
- 21.5" touch screen technology
- Superior HD image resolution for improved procedure accuracy, speed and physician confidence level
- Patented probe design
- Sagittal array provides for 140 mm length of view encompassing the bladder, prostate and perineal tissue
- Advanced modular software design provides for future upgrade path
- Convertible to desktop or portable without losing any functionality
- Dual battery backup built-in to protect from power loss

*Certain products shown are not available for sale currently.

Sonalis[®] Transducers

We offer probes for the most common types of procedures and imaging, all with the enhanced image quality and high resolution that you would expect from the **BEST** name in ultrasounds!

- 8L2A Linear Array 
- 12L5A Linear Array 
- 14L3 Linear Array 
- 15LW4 Linear Array 
- 15LA Linear Array 
- 15L4A Linear Array 
- 16L5 Linear Array 
- 8V3 Phased Array 
- 4V2A Phased Array 
- 5C2A Curved Array 
- 9MC3 Curved Array 
- 8EC4A Endocavity 
- XY-BI-Plane Phased Array 
- 10EC4A Endocavity 
- 10BP4 Bi-Plane 
- 8BP4 Bi-Plane 
- 8TE3 Trans-esophageal 
- Pedoff 
- 16HL7 High Frequency Linear Array 


TeamBest Global Companies ©2024

Best Nomos • 1 Best Drive • Pittsburgh, PA 15202 USA • tel 412 312 6700 or 800 70 NOMOS • fax 412 312 6701 • www.nomos.com

AccuBoost[®] DIGITAL

PARTIAL BREAST BRACHYTHERAPY SOLUTION

AccuBoost[®] Digital—A digital mammography-guided Non-Invasive Breast Brachytherapy (NIBB) equipment solution. This complete, turn-key solution provides all the components needed for this advanced treatment, bundled together in a streamlined offering.



HUESTIS MEDICAL
A TEAMBEST GLOBAL COMPANY

One Industrial Way | Tyngsboro, MA 01879

1-800-972-9222

www.accuboot.com
www.huestis.com
www.teambest.com



Welcome to Our World

TeamBest Global

www.teambest.com

*Ready to serve all your
needs globally*

NON-PROFIT

- **Best Cure Foundation**
www.bestcure.md
- **Brachytherapy Research & Educational Foundation**
www.brachytherapy.org
- **Global Best Cure Foundation**
www.globalbestcure.org

FOR-PROFIT

- **Arplay Medical**
www.arplay.com
- **Best ABT**
www.bestabt.com
- **Best Automation & Robotics**
www.teambestautomation.com

- **Best Cyclotron Systems, Inc.**
www.bestcyclotron.com
- **Best Dosimetry Services**
www.bestdosimetry.com
- **Best Entertainment**
www.bestentertainment.ca
- **Best Medical Canada**
www.bestmedical.ca or
www.mosfet.ca
- **Best Medical Capital, Inc.**
www.teambest.com
- **Best Medical International**
www.bestmedical.com
- **Best NOMOS**
www.nomos.com

- **Best Particle Therapy**
www.bestproton.com
- **Best Theratronics**
www.theratronics.com
- **Best Vascular (Novoste)**
www.bestvascular.com or
www.novoste.com
- **CNMC Company**
www.cnmcco.com
- **Huestis Medical**
www.huestis.com
- **Kitsault Energy**
www.kitsaultenergy.com

TEAMBEST GLOBAL®



www.teambest.com

Best Medical International
7643 Fullerton Road, Springfield, VA 22153
tel: 703 451 2378 800 336 4970

TeamBest Global Companies ©2024

TEAMBEST GLOBAL®



www.teambest.com